**SARS-CoV-2 Antigen Rapid Test Cassette**

Instruction for Use For Private Use  
(Self-testing)

**Package Information**

**REF CP01750011**

Components	1 T
Test Cassette	×1
Sample extraction solution	×1
Disposable sampling swab CE 0197 93/42/EEC	×1
Biohazard Specimen Bag	×1

**REF CP01760011**

Components	20 T
Test Cassette	×20
Sample extraction solution	×20
Disposable sampling swab CE 0197 93/42/EEC	×20
Biohazard Specimen Bag	×20

**REF CP01770011**

Components	50 T
Test Cassette	×50
Sample extraction solution	×50
Disposable sampling swab CE 0197 93/42/EEC	×50
Biohazard Specimen Bag	×50

**REF CP01800011**

Components	1 T
Test Cassette	×1
Sample extraction solution	×1
Specimen Collection Swab Model N150 CE REGULATION (EU) 2017/745	×1
Biohazard Specimen Bag	×1

**REF CP01810011**

Components	20 T
Test Cassette	×20
Sample extraction solution	×20
Specimen Collection Swab Model N150 CE REGULATION (EU) 2017/745	×20
Biohazard Specimen Bag	×20

**REF CP01820011**

Components	50 T
Test Cassette	×50
Sample extraction solution	×50
Specimen Collection Swab Model N150 CE REGULATION (EU) 2017/745	×50
Biohazard Specimen Bag	×50

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**Test Cassette**

Rapid Test for Detection of SARS-CoV-2's Nucleocapsid protein in white plastic cassette packed in aluminum foil bag. For single use only.

**Sample extraction solution**

300 µL per test.

**Others**

-Instructions for use

**Materials required but not provided**

Clock

**Intended Use**

The cassette is used for qualitative detection of antigen of SARS-CoV-2 in human nasal swabs.

The cassette is intended for screening of patients suspected for infection with SARS-CoV-2, and as an aid in the diagnosis of the coronavirus disease (COVID-19).

For single use only. The test cassette can be used from the age of 18 years. Children and adolescents under the age of 18 may only use the test cassette under adult supervision.

Our test detected the nucleocapsid protein of SARS-CoV-2 and not the surface protein ("spike"). Therefore, the test also recognizes the British variant.

For self-testing.

**Summary**

COVID-19 is a respiratory disease caused by infection with SARS-CoV-2 virus. Common signs of infection include respiratory symptoms, fever, cough, shortness of breathing difficulties. In severe cases, infection can cause pneumonia, severe acute respiratory syndrome (SARS), kidney failure and death.

**Method**

Gold Immunochromatographic Assay (GICA).

**Principle**

→ The Cassette is a one-step lateral flow chromatographic immunoassay. The test strip in the device includes: 1) a conjugate pad containing anti-SARS-CoV-2's Nucleocapsid protein antibody, Chicken IgY antibody, all of which are conjugated to colloidal gold, and 2) nitrocellulose membrane containing a test line (T line) and a control line (C line).

→ The T line is coated with anti-SARS-CoV-2's Nucleocapsid protein antibody, when the specimen is added, it migrates by capillary diffusion rehydrating the gold conjugate, if present in specimen, SARS-CoV-2's Nucleocapsid protein and its antibody labeled with colloidal gold formed antigen-antibody complexes. These complexes will continue to migrate along the strip until the T line, where they are captured by the anti-SARS-CoV-2's Nucleocapsid protein antibody generating a visible red violet line in T line. If the specimen does not contain SARS-CoV-2 or the SARS-CoV-2 level is below the lower level, the T line will not appear.

→ The C line is coated with Goat anti-chicken antibodies, which should bind to the gold-chicken IgY antibodies conjugate and form a red violet line regardless of the presence of SARS-CoV-2's Nucleocapsid protein.

**Storage and validity**

1. Sealed in aluminum foil bag at 2-30 °C .Valid for 24 months.
2. Protect from light and don't freeze.

**Specimen transport and storage**

Type: Nasal Swabs

1. The samples should be used as soon as possible after collected within 1 hour at room temperature (15-30°C).
2. After processing the swab in the extraction solution, the sample should be analyzed within 30 minutes.

**Test procedure**

Do not open the pouch until you are ready to perform a test, and the single-use test is suggested to be used under low environment humidity (RH≤70%) within 1 hour.


**Preparation**

1. Allow all cassette components and specimens to reach room temperature between 15°C~30°C prior to testing.
2. Remove the test cassette from the foil pouch and place on a clean dry surface.


**Operation**

1. Specimen collection: Carefully insert the swab into one nostril, the swab tip should be inserted up to 2.5 cm (1 inch) from the edge of the nostril. Roll the swab 5 times along the mucosa inside the nostril to ensure that both mucus and

cells are collected. Using the same swab, repeat this process for the other nostril.


2. Carefully tear off the sealing paper on the mouth of the sample extraction liquid tube, being careful not to spill the liquid.


3. Withdraw the swab from the nasal cavity, completely immerse the swab head in the sample extraction buffer in the tube.

4. Vigorously mix the solution by rotating the swab forcefully against the side of the tube at least 10 times (while submerged) for about 1 minute

5. Squeeze the tube 5 times by hand to ensure that the sample on the sampling swab is fully eluted into the sample extraction solution. (Tips: Squeeze the swab head along the inner wall of the extraction tube)

6. Cover the drip head to mix the liquid thoroughly

7. Dispense 3 drops (80µL) of the specimen into the circular sample well on the cassette. Interpret the test results at 15~20 minutes. (Tips: Do not interpret the results after 20 minutes.)


8. After the test is completed, place the test cassette / swab / solution in the Biohazard Specimen Bag and dispose all test kit materials via biohazard waste disposal protocol. Lastly wipe down all surfaces the kit came in contact with and wash your hands or use hand sanitizer

**Interpretation of results**

1. Positive: Both of T line and C line appear, as shown in the figure.
2. Negative: T line do not appear, while C line appear, as

shown in the figure.

3. Invalid: If C line does not appear, the result is invalid whether the T line appear or not, as shown in the figure. Prepare another cassette and specimen to test again. Tips: please analysis the volume of sample addition.

**What to do with a positive result?**

- There is currently a suspicion of COVID-19 infection
- Immediately contact a doctor/family doctor or the local health authority
- Comply with local self-insulation guidelines
- to have a PCR confirmation test carried out

**What should I do if the result is negative?**

The rules that apply locally to contain the pandemic with regard to contact with other people and other protective measures continue to apply.

A COVID-19 infection can also be present if the rapid test result is negative.

In case of suspicion, a rapid test should be repeated after 1 to 2 days, as an infection cannot be reliably detected in all phases of the infection using the rapid test.

**What should I do if the result is invalid?**

An invalid test result can be caused by an incorrect test execution.

You should repeat the test.

If invalid test results continue to occur despite strict adherence to the instructions for use, you should contact a doctor or a COVID-19 test center and discuss how to proceed.

**Benefit of the testing process**

Rapid antigen tests (also referred to as point-of-care tests, POCT) can be carried out outside of the laboratory near the patient and faster than the PCR (polymerase chain reaction) in 15-30 minutes. Suitable SARS-CoV-2 POC antigen tests can therefore play a role in situations in which a quick result is important and in which the contagiousness (infectiousness) of people should be assessed promptly and on site. The use of rapid antigen tests is so intended to interrupt the contribute to transmission through targeted isolation of the infected and their close contacts.

**Limitation of Test Methods**

1. The result of the product should not be taken as a confirmed diagnosis, for clinical reference only.

Judgment should be made along with RT-PCR results, clinical symptoms, epidemiological information and further clinical data.

2. The contents of this cassette are to be used for the qualitative detection of SARS-CoV-2 antigens from nasal swab.
3. Sensitivity of the test after the first seven days of the onset of symptoms has been demonstrated to decrease as compared to a RT-PCR SARS-CoV-2 assay.
4. The performance of this test has not been evaluated for use in patients without signs and symptoms of respiratory infection and performance may differ in asymptomatic individuals.
5. The Sample buffer and test cassette must be equilibrated to room temperature (15 °C~30 °C) before used, otherwise the results may be incorrect.
6. A negative test result may occur if the level of antigen in a sample is below the detection limit of the test or if the sample was collected improperly.
7. Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.
8. React less than 15 minutes may lead a false negative result; React more than 20 minutes may lead a false positive result.
9. Positive test results do not rule out co-infections with other pathogens.
10. Clinical performance was evaluated with frozen samples, and performance may be different with fresh samples.
11. Users should test specimens as quickly as possible after specimen collection.
12. After processing the swab in the extraction solution, the sample should be analyzed within 30 minutes.

**Performance Data**

**1. Clinical validation**

**Sample type: nasal swab**

The results have been summarized in the following tables:

Method	PCR		Total Results
	Positive	Negative	
SARS-Cov-2 Antigen Rapid Test	125	1	126
	7	128	135
<b>Total Results</b>	<b>132</b>	<b>129</b>	<b>261</b>

Positive coincidence rate =94.70%, (95% CI: 89.46%-97.41%);

Negative coincidence rate =99.22%, (95% CI: 95.74%-99.86%);

Total coincidence rate =96.93%, (95% CI: 94.07%-98.44%)

**2. Sensitivity**

The sensitivity was tested with three reference materials. These reference materials differed in their concentration of inactivated SARS-CoV-2 virus culture. The

concentrations of the reference material used were:

S1 = 0.8 × 10<sup>2</sup> TCID<sub>50</sub> / mL

S2 = 1.6 × 10<sup>2</sup> TCID<sub>50</sub> / mL

S3 = 3.2 × 10<sup>2</sup> TCID<sub>50</sub> / mL

The S1 reference material gave a negative test result.

The S2 and S3 reference material each produced a positive test result.

As a result, the minimum detection limit of the test cassette is 1.6 × 10<sup>2</sup> TCID<sub>50</sub> / mL

**3. Precision**

1 precise reference material from the manufacturer, consisting of inactivated SARS-CoV-2 virus culture and a defined concentration, was tested 10 times individually. The depth of color development of the test line (T) should be uniform and consistent, and the results should all be positive. A uniform, consistent color and a positive test result were achieved in each of the 10 tests.

The concentration of the reference material used was 4.8 × 10<sup>2</sup> TCID<sub>50</sub> / mL.

**4. Cross-reactivity**

Cross-reactivity of the cassette was evaluated. The results showed no cross reactivity with the following microorganism.

No.	Microorganism	Conc.	Cross reactivity (YES NO)
1	Human Coronavirus OC43	1×10 <sup>5</sup> pfu/mL	NO
2	Human Coronavirus 229E	1×10 <sup>5</sup> pfu/mL	NO
3	Human Coronavirus NL63	1×10 <sup>5</sup> pfu/mL	NO
4	Influenza A H1N1 (2009)	1×10 <sup>5</sup> pfu/mL	NO
5	MERS-coronavirus	1×10 <sup>5</sup> pfu/mL	NO
6	SARS-coronavirus	1×10 <sup>5</sup> pfu/mL	NO
7	Influenza A H3N2	1×10 <sup>5</sup> pfu/mL	NO
8	Influenza B Yamagata	1×10 <sup>5</sup> pfu/mL	NO
9	Influenza B Victoria	1×10 <sup>5</sup> pfu/mL	NO
10	Respiratory syncytial virus A	1×10 <sup>5</sup> pfu/mL	NO
11	Respiratory syncytial virus B	1×10 <sup>5</sup> pfu/mL	NO
12	Adenoviridae 1	1×10 <sup>5</sup> pfu/mL	NO
13	Adenoviridae 2	1×10 <sup>5</sup> pfu/mL	NO
14	Adenoviridae 3	1×10 <sup>5</sup> pfu/mL	NO
15	Adenoviridae 4	1×10 <sup>5</sup> pfu/mL	NO
16	Adenoviridae 5	1×10 <sup>5</sup> pfu/mL	NO
17	Adenoviridae 7	1×10 <sup>5</sup> pfu/mL	NO
18	Adenoviridae 55	1×10 <sup>5</sup> pfu/mL	NO
19	Enterovirus EV71	1×10 <sup>5</sup> pfu/mL	NO
20	Enterovirus CA16	1×10 <sup>5</sup> pfu/mL	NO
21	Enterovirus CA10	1×10 <sup>5</sup> pfu/mL	NO
22	Enterovirus CB5	1×10 <sup>5</sup> pfu/mL	NO
23	Enterovirus CA24	1×10 <sup>5</sup> pfu/mL	NO
24	Enterovirus CB4	1×10 <sup>5</sup> pfu/mL	NO
25	Enterovirus CB3	1×10 <sup>5</sup> pfu/mL	NO

26	Enterovirus CB2	1×10 <sup>5</sup> pfu/mL	NO
27	Enterovirus CB1	1×10 <sup>5</sup> pfu/mL	NO
28	Enterovirus CA6	1×10 <sup>5</sup> pfu/mL	NO
29	EB virus	1×10 <sup>5</sup> pfu/mL	NO
30	Human cytomegalovirus	1×10 <sup>5</sup> pfu/mL	NO
31	Human Rhinovirus A30	1×10 <sup>5</sup> pfu/mL	NO
32	Human Rhinovirus A31	1×10 <sup>5</sup> pfu/mL	NO
33	Human Rhinovirus A2	1×10 <sup>5</sup> pfu/mL	NO
34	Human Rhinovirus A81	1×10 <sup>5</sup> pfu/mL	NO
35	Human Rhinovirus B52	1×10 <sup>5</sup> pfu/mL	NO
36	Human Rhinovirus B70	1×10 <sup>5</sup> pfu/mL	NO
37	Human Rhinovirus B72	1×10 <sup>5</sup> pfu/mL	NO
38	Metapneumovirus A2	1×10 <sup>5</sup> pfu/mL	NO
39	Metapneumovirus Type B1	1×10 <sup>5</sup> pfu/mL	NO
40	Metapneumovirus Type B2	1×10 <sup>5</sup> pfu/mL	NO
41	Measles virus	1×10 <sup>5</sup> pfu/mL	NO
42	Rubella virus	1×10 <sup>5</sup> pfu/mL	NO
43	Mumps virus	1×10 <sup>5</sup> pfu/mL	NO
44	Boca virus	1×10 <sup>5</sup> pfu/mL	NO
45	Parainfluenza Virus 1	1×10 <sup>5</sup> pfu/mL	NO
46	Parainfluenza Virus 2	1×10 <sup>5</sup> pfu/mL	NO
47	Parainfluenza Virus 3	1×10 <sup>5</sup> pfu/mL	NO
48	Parainfluenza Virus 4	1×10 <sup>5</sup> pfu/mL	NO
49	Bordetella pertussis	1×10 <sup>6</sup> cfu/mL	NO
50	Candida albicans	1×10 <sup>6</sup> cfu/mL	NO
51	Legionella pneumophila	1×10 <sup>6</sup> cfu/mL	NO
52	Haemophilus influenzae	1×10 <sup>6</sup> cfu/mL	NO
53	Human Metapneumovirus	1×10 <sup>6</sup> cfu/mL	NO
54	Streptococcus pneumoniae	1×10 <sup>6</sup> cfu/mL	NO
55	Streptococcus pyogenes	1×10 <sup>6</sup> cfu/mL	NO
56	Mycobacterium tuberculosis	1×10 <sup>6</sup> cfu/mL	NO
57	Pneumocystis jirovecii (PJP)	1×10 <sup>6</sup> cfu/mL	NO
58	Staphylococcus aureus	1×10 <sup>6</sup> cfu/mL	NO
59	Staphylococcus epidermidis	1×10 <sup>6</sup> cfu/mL	NO
60	Mycoplasma pneumoniae	1×10 <sup>6</sup> cfu/mL	NO
61	Chlamydia pneumoniae	1×10 <sup>6</sup> cfu/mL	NO

**5. Interference Substances**

The test results do not be interfered with the substance at the following concentration:

No.	Substances	Conc.
1	Ibuprofen	1mg/mL
2	Tetracycline	3µg/mL
3	Chloramphenicol	3µg/mL
4	Erythromycin	3µg/mL
5	Tobramycin	5%
6	Throat spray (Menthol)	15%
7	Mupirocin	10mg/mL
8	Throat lozenge (Menthol)	1.5mg/mL
9	Oseltamivir	5mg/mL
10	Naphthoxoline hydrochloride nasal	15%

	drops	
11	Mucin	0.50%
12	Fisherman's Friend	1.5mg/mL
13	Compound Benzocain Gel	1.5mg/mL
14	Cromoglycate	15%
15	Phenylephrine Hydrochloride	15%
16	Afrin (Oxymetazoline)	15%
17	Fluticasone propionate spray	15%
18	Whole Blood	4%

**Precautions**

- For in vitro diagnostic use only.
- Do not use package when damaged, unclear label or even expired.
- The test shall be operated in strict accordance with the instructions.
- The results must be interpreted within 15 minutes.
- The disposable cassette is referred to biological waste when used. (e.g. in a tear-resistant, sealed plastic bag in the residual waste)
- Dispose the cassette and components according to local guidelines.
- Ensure that the test is performed hygienically so as not to contaminate the swab, reagents or test cassette. This might lead to incorrect test results.
- Make sure that the test cassette or other components do not get wet. Otherwise, this may lead to incorrect results.
- Avoid splashing specimen material. This can lead to contamination and, in the case of positive samples, to the risk of infection.
- Avoid splashing, ingestion, skin and mucous membrane and eye contact of/with buffer. This could lead to irritation.
- Wash your hands after using the extraction solution or rinse the affected parts of the body thoroughly with water if it comes into contact with the extraction solution. If symptoms persist, you should seek medical advice / assistance.
- Do not touch the test window or sample window of the test cassette to avoid contamination and false results.
- Do not mix buffer solutions or other components of different batches. This might lead to incorrect results.
- Do not mix patient samples. This might lead to incorrect results.

**References**

- 1.Chen Wei, Zhang Chunyang, Zhu Ying, Zhang Yanhua, You Libin, Wu Bingshan, Huang Zhimiao, Zheng Kuicheng, Weng Yuwei. Comparison of nucleic acid detection in pharyngeal swabs and sputum specimens in 4 cases of new coronavirus infection. Chinese Journal of Zoonoses 2020; 1-7 (in Chinese)
- 2.Guo Yuanyuan, Wang Kun, Zhang Yu, Zhang Wenjia,

Wang Liying, Liao Yi. Comparison and analysis of detection performance of six domestic new coronavirus nucleic acid detection reagents. Chongqing Yixue 2020; 1-10 (in Chinese).

**Merlin Biomedical (Xiamen) Co., Ltd.**

Add: 4th Floor, Building B3,  
2054 West Wengjiao Rd.,  
Haicang, Xiamen, 361028, China  
Tel: 86-592-5210772; 86-592-5210773  
Fax: 86-592-5210772  
Email: info@merlinbio.com.cn  
Website: http://www.merlinbio.com.cn

**Qarad EC-REP BV**

Pas 257, 2440 Geel, Belgium

**Distributor & Importer:**

**Arterial Global Sdn. Bhd. (894261-H)**

3-2, Jalan Ara SD 7/3B,  
Bandar Sri Damansara,  
52200 Kuala Lumpur, Malaysia  
Tel/Fax: 03-6263-1491  
Email: sales@arterialglobal.com.my

**Instruction Approval and Revision Date**

Approval Date: 12 Aug 2020 en-Self-testing

Revision Date: 02 Mar 2022 en-Self-testing


Date of Issue: 02 Mar 2022 en-Self-testing

**Symbols**

	Catalogue number
	Consult instructions for use
	In vitro diagnostic medical device
	Manufacturer
	Temperature limit
	Batch code
	Use-by date
	Keep dry
	Authorized representative in the European Community
	Don't use the product when the package is damaged
	Do not re-use
	Warning
	Contains sufficient for <n> tests
	Date of manufacture
	Sterilized using ethylene oxide
	Do not re-sterilize

	Avoid overexposure to the sun
	Test Cassette
	Sample extraction solution
	Disposable sampling swab
	Specimen Collection Swab
	Biohazard Specimen Bag

**Reporting COVID-19 Result Method**

**QR code**

Please scan the QR code on the box for more details

**Disposal Procedure**

Throw all Test Cassette, Sample extraction solution, Tube, Dropper, swab into Biohazard Specimen Bag. Throw the Biohazard Specimen Bag and others general waste into general waste bin.

Note: All components of the kit: Discard as per the Guidelines for handling, Treatment and Disposal of waste Generated during Treatment/Diagnosis/Quarantine of COVID-19 patients issued by Local Regulatory Requirements or Regulatory Authorities.

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11	Mucin	0.50%
12	Fisherman's Friend	1.5mg/mL
13	Compound Benzocain Gel	1.5mg/mL
14	Cromoglycate	15%
15	Phenylephrine Hydrochloride	15%
16	Afrin (Oxymetazoline)	15%
17	Fluticasone propionate spray	15%
18	Whole Blood	4%

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### Rujukan

- Chen Wei, Zhang Chunyang, Zhu Ying, Zhang Yanhua, You Libin, Wu Bingshan, Huang Zhimiao, Zheng Kuicheng, Weng Yuwei. Comparison of nucleic acid detection in pharyngeal swabs and sputum specimens in 4 cases of new coronavirus infection. Chinese Journal of Zoonoses 2020; 1-7 (in Chinese)
- Guo Yuanyuan, Wang Kun, Zhang Yu, Zhang Wenjia, Wang Liying, Liao Yi. Comparison and analysis of detection

performance of six domestic new coronavirus nucleic acid detection reagents. Chongqing Yixue 2020; 1-10 (in Chinese).

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
### Melaporkan Kaedah Keputusan COVID-19

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